

## Fresh Product

### Whole Peripheral Blood Vacutainer<sup>®</sup> Tube Collection, Hodgkin's Lymphoma

Catalog#	PB010F-HL	10 mL
	PB020F-HL	20 mL
	PB030F-HL	30 mL
	PB040F-HL	40 mL
	PB050F-HL	50 mL
	PB060F-HL	60 mL
	PB070F-HL	70 mL
	PB080F-HL	80 mL
	PB090F-HL	90 mL
	PB100F-HL	100 mL

## Product Description

Human Whole Peripheral Blood is collected from Hodgkin's lymphoma patients in 10 mL sterile vacutainer(s) containing EDTA.

Fresh products have a high viability without the detrimental effects of freezing, thawing, and exposure to cryoprotectants.

Blood is obtained using Institutional Review Board (IRB) approved consent forms and protocols.

## Sample Collection and Processing

All samples are collected on-site at our Stem Cell Collection Center. Samples are quickly processed in our on-site laboratory to achieve maximum viability and quality.

## Format

Fresh whole peripheral blood is shipped at ambient temperature.

Specific containers and media can also be prepared as requested by the customer.

## Storage

Fresh products should be used or processed immediately upon receipt. The warranty only covers items whose specifications are tested at the time they are received.

## Warning

This product contains human tissue or other biological material and MUST be handled at Biosafety Level 2 or higher. All biological products should be treated as potentially infectious or contaminated material, even if infectious disease screening reports are negative. Follow universal precautions and wear appropriate personal protective equipment.

## Product Warranty

For our product warranty, please review our Terms and Conditions at [stemexpress.com/terms-and-conditions/](http://stemexpress.com/terms-and-conditions/).

FOR IN VITRO RESEARCH USE ONLY. NOT APPROVED FOR HUMAN OR VETERINARY USE IN VIVO.